



2024/2025 Catalogue



*Sutures Designed for
All Animal Surgeries*

Sutures *Non-Absorbable*

SHORTY™ Sutures - Equine and Small Animal

PV		POLYVINYLIDENE FLUORIDE, NON-ABSORBABLE						
ABSORPTION PROFILE			Exceptional Knot Security, knots flatten and lock against themselves, unwetted by blood, inert Non-Absorbable					
ORDER CODE	USP GAUGE	EP	NEEDLE TYPE	NEEDLE SIZE	NEEDLE	SUTURE LENGTH	QTY/BOX	PRICE
PV-20425D	2-0	3	3/8 Reverse Cutting ▽	25mm		45cm	36	
PV-30425D	3-0	2	3/8 Reverse Cutting ▽	25mm		45cm	36	


Double forward throw

(Surgeons Throw)



Single reverse throw

Creating a square knot, with which to adjust tension.



Single forward throw

(Locking the knot)



Final result

9 interrupted sutures, with long tails, from a 45cm length.



Attention All Veterinarians . . .



SHORTY™ Sutures released!

Ever just needed a short length high vis suture for a small incision?




Look no further than **SHORTY™**.

- Low Cost - No Waste.



Innovation: Vets Talk, SilverGlide Listens






Internal Sutures

PDO		Polydioxanone Sutures, Purple Monofilament, Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE	
	PO-CLK-01430D	1	50mm	3/8 Circle Reverse Cutting	45cm	36		
	PO-01130D	1	50mm	3/8 Circle Reverse Cutting	100cm	36		
	PO-2075D	2-0	26.5mm	1/2 Circle Taper Point	70cm	36		
	PO-2078D	2-0	36mm	1/2 Circle Taper Point	70cm	36		
	PO-0018D	0	36mm	1/2 Circle Taper Point	100cm	36		

Skin Closure



N		Nylon Sutures, Black Monofilament, Non-Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY/ BOX	PRICE	
	N-20710D	2-0	30mm	3/8 Circle Reverse Cutting	70cm	36		
	N-00110D	0	30mm	3/8 Circle Reverse Cutting	100cm	36		

Skin Closure







SMF		Supramid, Fluro Yellow, Non-Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY/ BOX	PRICE	
	SMF-3077D	3-0	24mm	3/8 Circle Reverse Cutting	70cm	36		
	SMF-00113D	0	37mm	1/2 Circle Reverse Cutting	100cm	36		
	SMF-01130D	1	50mm	3/8 Circle Reverse Cutting	100cm	36		

Equine Surgery

* Low Cost Sutures for Small Incisions



RR		Polypropylene Fluro Pink Monofilament, Non-Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY/ BOX		
	RR-CLK-00430D	0	50mm	3/8 Circle Reverse Cutting	45cm	36		
	RR-CLK-02430D	2	50mm	3/8 Circle Reverse Cutting	45cm	36		
	RR-02130D	2	50mm	3/8 Circle Reverse Cutting	100cm	36		

Excellent contrast for suture visualisation during suture placement and removal

FG		ORANGE POLYVINYLIDENE FLUORIDE , NON-ABSORBABLE						
NEEDLE	ORDER CODE	USP	NEEDLE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE	
	FG-3077D	3-0	24mm	3/8 Circle Reverse Cutting	70cm	36		
	FG-20760D	2-0	60mm	Straight Cutting	70cm	36		
	FG-01930D	4	50mm	3/8 Circle Reverse Cutting	90cm	24		
	FG-20710D	2-0	30mm	3/8 Circle Reverse Cutting	70cm	36		
	FG-00913D	3-5	37mm	1/2 Circle Reverse Cutting	90cm	24		
	FG-00930D	3-5	50mm	3/8 Circle Reverse Cutting	90cm	24		

* Race track dispenser pack

Internal Sutures

PGA		Polyglycolic Acid Sutures, Braided Violet, Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY		
	PG-0018D	0	36mm	1/2 Circle Taper Point	100cm	36		
	PG-CLK-02430D	2	50mm	3/8 Circle Reverse Cutting	45cm	36		

Examples of high visibility sutures making for easier suture placement, location and removal.




Fetlock Repair




Arthroscopy Port Closure






Skin Closure

N		Nylon Sutures, Black Monofilament, Non-Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY/ BOX	PRICE	
	N-01130D	1	50mm	3/8 Circle Reverse Cutting	100cm	36		

SMF		Supramid, More Fluro Yellow, Non-Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY/ BOX	PRICE	
	SMF-031360N	3	65mm	Straight Cutting	130cm	12		
	SMF-031375N	3	75mm	1/2 Circle Conventional Cutting Straight Shank	130cm	12		
	SMF20710D	2.0	30mm	3/8 Circle Reverse Cutting	70cm	36		
	SMF00110D	0	30mm	3/8 Circle Reverse Cutting	100cm	36		

Internal Closure

GCL		Glycolide & E-Cap Sutures, Purple Monofilament, Absorbable						
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE	
	EC-011339N *	1	65mm	1/2 Circle Taper Point	130cm	12		

* Bovine Uterus closure



PGA Polyglycolic Acid Sutures, Braided Violet, Absorbable							
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE
	PG- 3075D	3-0	26.5mm	1/2 Circle Taper Point	70cm	36	
	PG- 2075D	2-0	26.5mm	1/2 Circle Taper Point	70cm	36	
	PG -021339N	3	65mm	1/2 Circle Taper Point	130cm	12	


Skin Closure

PDO Polydioxanone Sutures, Purple Monofilament, Absorbable							
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE
	PO-4074D	4-0	22mm	1/2 Circle Taper Point	70cm	36	
	PO-3075D	3-0	26.5mm	1/2 Circle Taper Point	70cm	36	
	PO-0018D	0	36mm	1/2 Circle Taper Point	100cm	36	
	PO-021339N*	2	65mm	1/2 Circle Taper Point	130cm	12	
	PO-3077D	3-0	24mm	3/8 Circle Reverse Cutting	70cm	36	
	PO-20726D	2-0	26mm	3/8 Circle Reverse Cutting	70cm	36	




Needle Shapes

1/2 Circle




3/8 Circle




Straight Cutting



1/2 Circle Conventional Cutting /
Straight Shank (Ski)




Needle Types




Round bodied taper point needle ◉



Reverse cutting needle ▽












Conventional cutting needle △



Small Animal Surgery

Internal Closure -







GCL Glycolide & E-Cap Sutures, Purple Monofilament, Absorbable							
NEEDLE	CODE	USP GAUGE	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY	PRICE
	EC-3075D	3-0	26.5mm	1/2 Circle Taper Point	70cm	36	
	EC-2075D	2-0	26.5mm	1/2 Circle Taper Point	70cm	36	
	EC-50416N	5-0	16mm	3/8 Circle Reverse Cutting	45cm	12	
	EC-4076D	4-0	18mm	3/8 Circle Reverse Cutting	70cm	36	
	EC-3077D	3-0	24mm	3/8 Circle Reverse Cutting	70cm	36	
	EC-20726D	2-0	26mm	3/8 Circle Reverse Cutting	70cm	36	




PGF Polyglycolic Acid Fast Absorbing Sutures - Dental							
NEEDLE	USP GAUGE	CODE	EP	NEEDLE SIZE	NEEDLE TYPE	SUTURE LENGTH	QTY
	5-0	PGF-50413N	1	13mm	1/2 Circle Reverse Cutting	45cm	12
	4-0	PGF-40422N	1.5	16mm	3/8 Circle Reverse Cutting	45cm	12
	3-0	PGF-3046N	2	18mm	3/8 Circle Reverse Cutting	45cm	12





Suture Material Absorption Profile

NON-ABSORBABLE		Material Description	Qualities	SIMILAR TO
NYLON		Black Monofilament	High strength retention	Ethilon™ Monosof™ ^{1/2} Dermalon™ ¹ Dafilon™ ¹ Riverlon
SUPRAMID		Yellow polyamide multifilament	Easy handling, Ultra high visibility Fluorescent Yellow	
SUPRAMID		White polyamide multifilament	Easy handling	
BLUE POLYPROPYLENE		Dark blue monofilament	Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security	Prolene™ Surgipro™ ^{1/2} /Surgilene™ ^{1/2} Premilene™ ¹
FLUORESCENT PINK POLYPROPYLENE		Pink monofilament	Ultra high visibility, Fluorescent Pink Unwetted by blood, inert, knots flatten and lock against themselves for improved knot security.	Visipro™
POLYESTER		Braided polyester green silicone coated	High strength braid	Ticon™ Ethibond™

ABSORBABLE		Material Description	BSR Profile	Mass Absorption	SIMILAR TO
GLYCOLIDE & E-CAP		Monofilament purple	In vivo strength retention 60% at 7 days 25% at 14 days	Essentially complete between 2-4 months	Monocryl™ Caprosyn™ ^{1/2} /Biosyn™ ¹ Monosyn™ ¹ MonoQ™
POLYGLYCOLIC ACID		Braided violet	In vivo strength retention 65% at 14 days 20% at 21 days	Essentially complete between 3-6 months	Vicryl™ Safil™ ¹ Polisorb™ Villet™
POLYDIOXANONE		Monofilament purple	In vivo strength retention 80-90% at 14 days 60-80% at 28 days 50-80% at 42 days	Essentially complete between 6-9 months	PDS11™ MonoPlus™ ¹ Maxon™ PDM™

